




DATA SHEET

ENDOSCOPE WASHER MOD.

GANDY-90





SUMMARY

| | |
|--|-----------|
| 1. Introduction | 3 |
| 1.1 Preliminary Information | 3 |
| 1.2. Certifications and Standards | 3 |
| 1.3 Field of Application | 4 |
| 2. Device features | 4 |
| 2.1 Materials | 4 |
| 2.2 Dimensional parameters | 5 |
| 2.3 Structural features | 6 |
| 2.4 Functional features | 7 |
| 2.5 Tracking | 8 |
| 2.6 Reconditioning programs | 9 |
| 2.6.1 Time and consumption of individual phases | 9 |
| 2.6.2 Washing and disinfection features | 10 |
| 2.7 Safety | 11 |
| 2.7.1 General safety measures | 11 |
| 2.7.2 Alarms | 13 |
| 3. Installation requirements | 13 |
| 3.1. Electrical requirements | 13 |
| 3.2 Internet connection requirements | 14 |
| 3.3. Water requirements | 14 |
| 3.4 Environmental requirements | 14 |
| 3.5 Installation scheme | 15 |
| 4. Consumables and accessories | 15 |
| 4.1 Consumables | 15 |
| 4.2 Accessories | 16 |

1. Introduction

1.1 Preliminary Information

| MODEL | GANDY-90 |
|--------------------|---|
| Producer | NUOVA SB SYSTEMA S.R.L. |
| Manufacturer | Nuova SB System S.r.l. |
| CND | Z12011311 – STERILISERS FOR ENDOSCOPES |
| RDM | 2168938 |
| Designation of use | Washing, high-level disinfection and low-temperature chemical sterilization of thermolabile endoscopes. |
| Risk class | Class IIb |

1.2. Certifications and Standards

| | |
|----------------------------------|--|
| CE CERTIFICATE | CE Certificate - MED 31114 extended in the transitional period until the transition to the new medical device regulation MDR 2017/745 as specified in Art. 120.3c of the MDR as amended by Regulation (EU) 2023/607. |
| REFERENCE STANDARDS | |
| UNI EN ISO 15883-1: 2014 | Washer-disinfectors - Part 1: General requirements, terms and definitions and tests |
| UNI EN ISO 15883-4: 2019 | Washer-disinfectors - Part 4: Requirements and tests for washer-disinfectors using chemical disinfection for thermolabile endoscopes |
| UNI EN ISO 15883-5: 2021 | Washer-disinfectors - Part 5: Performance requirements and criteria for test methods to demonstrate cleaning effectiveness |
| UNI EN ISO 14937 | Sterilization of medical products - General requirements for the characterization of a sterilizing agent and for the development, validation and systematic control of a sterilization process for medical devices. |
| IEC 61010-1:2010 | Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements. |
| IEC 61010-2-040: 2020 | Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-040: Particular requirements for sterilizers and washing machines used to process medical materials. |
| EN 61326-1: 2022 | Electrical equipment for measurement, control and laboratory use - Electromagnetic compatibility requirements - Part 1: General requirements. |
| ISO 22196:2011 | Measurement of antibacterial activity on plastics and other non-porous surfaces. |
| UNI EN ISO 13850:2015 | Safety of machinery - Emergency stop function - design principles. |
| UNI CEI EN ISO 14971:2022 | Medical devices - Application of risk management to medical devices. |
| COMPANY STANDARDS | |
| UNI EN ISO 9001:2015 | Quality Management Systems - Requirements. |
| UNI CEI EN ISO 13485:2021 | Medical devices - Quality management systems - Requirements for regulatory purposes. |

1.3 Field of Application

The endoscope washing machine mod. GANDY-90 features two asynchronous washing tanks. The following table shows the types and number of endoscopes that can be reprocessed in each tank.

| ENDOSCOPE TYPE | NUMBER OF REPROCESSABLE ENDOSCOPES PER SINGLE CYCLE |
|--------------------------------|---|
| Flexible Endoscope | 1 |
| Bronchoscope | 2 |
| ENT without operative channels | 4 |
| ENT with one operative channel | 2 |

2. Device features

2.1 Materials

The **materials** conceived and used for the design and workmanship of the mod. GANDY-90 washer-disinfectors are shown in the following table.


| COMPONENT | MATERIAL |
|-------------------|--|
| Metal structure | Stainless Steel 316 |
| Fittings | Stainless Steel 316 |
| Tank | PMMA |
| Machine structure | Polyurethane coated with antibacterial paint |
| Internal tubes | Silicone or polyurethane |

This other table lists the general parameters of the equipment.


| PARAMETER | VALUE |
|---------------------------------|---------|
| Average noise emission | 60 db |
| Weight | 234 Kg |
| Maximum power input | 3.45 kW |
| Electrical safety class | I |
| Degree of electrical protection | IP45 |
| Type of biomedical device | B |

2.2 Dimensional parameters


Below there are the **views and relevant dimensional dimensions in millimeters** of the medical device mod. GANDY-90.




2.3 Structural features



Simple and intuitive 9.7" touch-screen software interface




Different coloured connectors and tubes for each channel for easy identification



2.4 Functional features

The following table shows the **functional characteristics** guaranteed by the machine during the endoscope reprocessing operating cycles for each tank. The various programs, with their durations and water consumption for each phase, will be described later in section 2.5.

| GUARANTEED FUNCTIONALITY DURING WASHING CYCLES | |
|--|--|
| Leak test performed at the beginning of the cycle and throughout the duration of the cycle. | |
| Automatic single-shot chemical dosing. | |
| Individual channel monitoring throughout the cycle for obstructions or disconnections. | |
| Pulsatile, unidirectional and self-regulating flow within each endoscope channel. | |
| Possibility to cycle in 'fast' mode by framing the endoscope and operator barcode. | |
| Possibility of programming self-disinfection cycles. | |
| Integrated capacitive system that heats the water to a suitable temperature to ensure perfect hot water washing. | |
| Possibility of discharging waste products into the sewage system. | |
| Immersion and spray washing and disinfection process. | |
| Automatic operator and endoscope recognition system via barcode gun. | |
| Real-time monitoring on the large touch-screen of all detected parameters (see figure below). | |



The other table shows the **operating functions** that the GANDY-90 washer-disinfector is capable of performing in addition to those normally guaranteed during endoscope instrument reprocessing cycles.

| EXTRA OPERATIONAL FEATURES |
|--|
| Three levels of accessibility: user, maintainer and admin. |
| Possibility of microbiological sampling directly from the wash tank. |
| Possibility to check the number of consumables in stock, with a warning message if stock is running low. |
| Monitoring the amount of chemical agent in the canisters. |
| Warning messages for filter change and routine maintenance. |

2.5 Tracking

The following table shows the **tracking** characteristics of the mod. GANDY-90 washer.

| TRACKING FEATURES |
|---|
| Complete tracking and recording of machine cycle data |
| Guaranteed paper tracking thanks to on-board thermal printer |
| The tracking of canisters and chemical agents used |
| Possibility of identifying whether a new canister you want to enter has already been used |
| Possibility of tracking thermal paper rolls |
| Operators, endoscopes and chemical agent canisters will be identified by their own barcode, which can be registered and identified on the equipment thanks to the barcode gun |
| Equipped with barcode gun. |
| Cycle traceability data can be exported via USB and Ethernet in .csv format. |

The **information and data** in the reports printed on thermal paper for each cycle to ensure paper tracking are shown in the following table.

| INFORMATION IN PAPER REPORTS |
|---|
| Hospital name. |
| Machine serial number. |
| Cycle number. |
| Selected washing program. |
| Date, time of loading, unloading of the instrument and total duration of the cycle. |
| Code of the responsible operator. |
| Endoscope serial number and its group. |
| Quantity of chemicals used and batch of canisters. |
| Phases of the cycle with relative duration. |
| Water circuit pressure and leak test. |
| Confirmation of cycle validation. |
| Possible red or orange alarm code with message of non-validation of the cycle. |
| Type of channels monitored. |

2.6 Reconditioning programs

The following section will specifically describe the washing and disinfection features within each reprocessing program, as well as detailing the timing of each step and the water consumption required in each step.

2.6.1 Time and consumption of individual phases

High disinfection programs

| | Program 1: Standard | Program 2: Standard double washing | Program 3: Advanced | Program 4: Disinfection phase | Program 6: Standard ENT | Program 9: Leak test repair |
|----------------------------|--------------------------------|---|--------------------------------|-------------------------------------|--------------------------------|-----------------------------------|
| Leak test | 20 s | 20 s | 20 s | 20 s | 20 s | 20 s |
| Pre-washing | 77 s | 77 s | 77 s | | 77 s | 77 s |
| Washing 1 | 355 s Contact time:300 s | 355 s Contact time:300 s | 355 s Contact time:300 s | | 355 s Contact time:300 s | 355 s Contact time: 300 s |
| Washing 2 | | 355 s Contact time:300 s | 355 s Contact time:300 s | | | |
| Rinse 1 | 95 s | 75 s | 95 s | | 95 s | 95 s |
| Rinse 2 | | 85 s | | | | |
| Disinfection | 355 s Contact time:300 s | 355 s Contact time:300 s | 505 s Contact time:450 s | 355 s Contact time:300 s | 355 s Contact time:300 s | 355 s Contact time:300 s |
| Final Rinse 1 | 65 s | 65 s | 75 s | 75 s | 65 s | 65 s |
| Final Rinse 2 | 65 s | 65 s | 65 s | 65 s | 65 s | 65 s |
| Final Rinse 3 | 85 s | 85 s | 85 s | 120 s | 85 s | 85 s |
| Drying | 87 s | 87 s | 87 s | 87 s | | 87 s |
| Total Time | 20 minutes | 27 minutes | 28 minutes | 13 minutes | 16 minutes | 19 minutes |
| Total water consumption | 46 L | 57 L | 50 L | 28 L | 46 L | 46 L |

In addition, the high disinfection program 11 is the one that allows for a fully customizable cycle.

Sterilization program

| | Program 1: Standard |
|-------------|---------------------------------|
| Leak test | 20 s |
| Pre-washing | 77 s |
| Washing 1 | 355 s Contact time: 300 s |
| Rinse 1 | 75 s |



| | |
|-------------------------|---------------------------------|
| Rinse 2 | 85 s |
| Disinfection | 655 s Contact time: 600 s |
| Final Rinse 1 | 65 s |
| Final Rinse 2 | 65 s |
| Final Rinse 3 | 85 s |
| Drying | 87 s |
| Total time | 26 minutes |
| Total water consumption | 53 L |

Other programs

| | Program 5: Drying | Program 7: Chemical self- disinfection | Program 8: Thermal self- disinfection | Program 10: Leak test |
|-------------------------|----------------------|--|---|-----------------------------|
| Leak test | | | | 120 s |
| Pre-washing | | 175 s | | |
| Disinfection | | 655 s Contact Time: 600 s | 2230 s Temperature from 65 to 82°C Contact time: 2160 s | |
| Final Rinse 1 | | 355 s | 115 s | |
| Final Rinse 2 | | 175 s | 125 s | |
| Final Rinse 3 | | 175 s | | |
| Final Rinse 4 | | 195 s | | |
| Drying | 200 s | | | |
| Total time | 200 s | 28 minutes | 41 minutes | 2 minutes |
| Total water consumption | | 42 L | 21 L | |

2.6.2 Washing and disinfection features

The following table lists the **chemical agents** manufactured and tested to ensure perfect disinfection of endoscopes in combination with the mod. GANDY-90 endoscope washer, including the **operating temperatures and mode of use**.



| WASHING | |
|---------------------------|--|
| Chemical used | Eco Zyme, multienzymes cleanser and class I medical device |
| Washing phase temperature | Between 30 and 45°C |
| Modalities of use | 1% dilution, 40 mL is dosed in 4 litres of water |
| DISINFECTION | |
| Chemical used | Ecosteril-F Plus, peracetic acid-based disinfectant and class IIb medical device |
| Washing phase temperature | <30°C, or at least the temperature of cold mains water |
| Modalities of use | 1% dilution, 70 mL is dosed in 7 litres of water |

2.7 Safety

The model GANDY 90 medical endoscope washer is equipped with all the necessary safety measures to guarantee the safety of the operator, endoscope, patient and machine itself, as well as to limit cross-contamination as much as possible. The following tables show the design solutions adopted to guarantee the safety of all those involved in the reprocessing cycle of flexible endoscopes.

2.7.1 General safety measures

Operator safety

| OPERATOR SAFETY SYSTEMS |
|--|
| Anti-crushing sensor active when closing the doors. |
| Closed-loop system to prevent leakage of fumes and chemicals into the environment. |
| Presence of an active carbon filter for each tank to remove toxic fumes from the use of peracetic acid disinfectant. |
| Locking of the doors during the entire cycle, which can only be opened at the end of the reconditioning cycle. |
| Doors with watertight closure. |
| Wash tank vent filter with 0.01 µm porosity |
| Caps on the detergent and disinfectant canisters integrated in the suction lances to facilitate changing canisters if the chemical agent is exhausted. |
| Different colored canister caps to distinguish detergent and disinfectant. |
| Emergency discharge in case of alarm and non-validated cycle. |
| Emergency unloading also possible from the main menu. |
| Possibility of safety rinsing of the tank. |

Safety for the endoscope

| ENDOSCOPE SAFETY SYSTEMS |
|--|
| Watertight doors closure with soft silicone sleeve. |
| Temperature monitored throughout the reconditioning cycle by temperature sensor. |
| Maximum permitted temperature of 45°C |
| Acoustic and visual alarm system capable of signaling possible channel disconnections or obstructions. |
| Continuous monitoring and measurement of all pressures generated in each individual endoscope channel. |

Safety for the machine

Safety systems dedicated to safeguarding the equipment, which then inevitably also affect endoscopes and operators, can be divided into **anti-corrosion and anti-contamination systems**.

| MACHINE SAFETY SYSTEMS |
|---|
| ANTI-CORROSION SYSTEMS |
| Materials resistant to corrosion caused by peracetic acid disinfectant (See section 2.1). |
| Scheduled periodic replacement of the most delicate components subject to the corrosive action of the peracetic acid disinfectant. |
| Cold disinfection to limit the corrosive effect of the peracetic acid disinfectant. |
| ANTI-CONTAMINATION SYSTEMS |
| Chemical and thermal self-disinfection. |
| Presence of non-return valves. |
| Disconnected water loading and unloading system. |
| Smooth equipment surfaces that can be easily washed with any commercially available detergent. |
| Antimicrobial coating of the entire instrumentation, endoscope washer and barcode gun, containing bacteriostatic additive in accordance with ISO 22196. |

General security systems





| GENERAL SAFETY SYSTEMS |
|--|
| Emergency button in accordance with ISO 13850. |
| Presence of IOlink sensors for monitoring machine parameters and working conditions. |
| Possibility of remote technical support. |

Limitation of cross-contamination

| CONTAMINATION LIMITATION SYSTEMS |
|---|
| CROSS-CONTAMINATION SYSTEMS |
| Hands-free technology: presence of an infrared sensor that detects the operator's foot, allowing the door to be opened without contact with the device. |
| Equipped with a stainless-steel basket compatible with the disinfection processes of the GANDY 90 washer-disinfector. |
| CYCLE ANTI-CONTAMINATION SYSTEMS |
| Inlet water filtration with double membrane filter with 0.2+0.1 µm porosity. |
| Air filtration used for drying with 0.01 µm filter. |

2.7.2 Alarms

The mod. GANDY-90 endoscope washer is equipped with **acoustic and visual** alarms capable of signaling to the operator any type of malfunction that may occur during the endoscope reprocessing cycle. The alarms are colour-coded according to the severity of the malfunction that occurs during the cycle.

| WASHING | |
|---|---|
| ALARM COLOUR CODE | SIGNIFICANCE |
|  | Visual warning only with message on the main page. |
|  | Acoustic and visual alarm, with message, which does not interrupt the cycle. |
|  | Acoustic and visual alarm, with message and error code, interrupting the cycle. |
|  | Acoustic and visual alarm, with message and error code, interrupting the cycle. |

3. Installation requirements

3.1. Electrical requirements

| ELECTRICAL REQUIREMENTS | |
|-------------------------|-------------|
| Supply voltage | 400V 3P + N |
| Power supply frequency | 50 Hz |
| Input power | 3.45 kW |
| Maximum amperage | 10 A |

3.2 Internet connection requirements

| INTERNET CONNECTION REQUIREMENTS | |
|----------------------------------|-----------------|
| LAN socket | LAN socket RJ45 |
| LAN socket data | IP address |
| | Subnet Mask |
| | DNS Address |
| | Gateway |

3.3. Water requirements


| REQUIREMENTS FOR WATER CONNECTION | |
|-----------------------------------|--|
| Water connections | Cold water inlet |
| | Drainage |
| WATER CHARACTERISTICS | |
| Minimum flow rate | 15 L/min |
| Hardness | < 120 mg/L of calcium carbonate (12°F) N.B.: In case of high-water hardness, with values above 12°F, the installation of a water softener is mandatory. For values close to, but below, 12°F, a water softener is strongly recommended. |
| Pressure | 2-4 bar |
| Temperature | <40°C If the inlet water temperature is between 30 and 40°C, the appliance will not use the internal boiler to heat the water further. |
| pH | 6.5-9 |
| Conductivity | At least 10 µS |

3.4 Environmental requirements

| ENVIRONMENTAL REQUIREMENTS | |
|--|-------------------|
| Room temperature | 5-60°C |
| Relative humidity | <80% |
| Maximum horizontal inclination allowed | 1° |
| Recommended ventilation frequency | > 10 Air change/h |

3.5 Installation scheme

INSTALLATION SHEET ENDOSCOPES-WASHER GANDY-90



4. Consumables and accessories

4.1 Consumables

| CONSUMABLE | CHARACTERISTICS |
|---------------------------|---|
| Ecosteril-F Plus | Peracetic acid disinfectant stored in 5-litre HDPE canister. |
| Eco Zyme | Multi-enzymatic cleaner stored in 5-litre HDPE canister. |
| Thermal paper KT 55 FA | Thermal paper with guaranteed 10-year stability. |
| Prefilter WFC-DPSHF4545 | Water prefilter with 0.45 µm porosity. |
| Prefilter WFCC2-DPPD10100 | Water prefilter with 1 µm porosity. |
| Filter WF4C2-DPSHF-2210 | Water filter with double filter membrane and 0.2+0.1 µm porosity. |
| Filter WSF-GPFL-0001-05 | Air filter with 0.01 µm porosity. |
| Filter M1000-8G-W | Leak test air filter with 0.01 µm filtration degree. |
| Filter 90DHV-PF0001EP | Tank vent air filter with 0.01 µm porosity. |
| Filter CamCarb CG ABS | Activated carbon filter. |

4.2 Accessories

| ACCESSORIES | CHARACTERISTICS |
|---------------------------|--|
| Connectors for endoscopes | Connectors are available for each channel of endoscopes of any model and brand, such as Pentax, Olympus, Fujinon, Storz, Wolf and others. |
| Stainless steel basket | Designed and manufactured to be contained within the wash tank, this basket can be used for handling endoscopes by limiting contact with them. |