



## THE C-EYE DEVICE

BY EMAGINE

CE marked for the following indications: keratoconus, ectasia after LASIK/PRK, pellucid marginal degeneration, infectious keratitis, sterile corneal melting, bullous keratopathy.



### Slit Lamp Mount

For modern CXL: fits on Haag Streit and Zeiss design-slit lamps.



### Range of Intensities

Delivers 3, 9, 15, 18 and 30 mW/cm<sup>2</sup>.



### Table Mount

For classic CXL: table mount accessory available.



### Battery-Operated

Charge via USB-C for maximal mobility.



### Thickness-Adjusted Profile

Beam profile delivers more energy in the corneal periphery.



### Pulsed and Continuous Light

For various keratoconus & infectious keratitis settings.



### Preset Protocols

Covering the most widely used CXL settings.

[www.emagine-eye.com](http://www.emagine-eye.com)

The C-Eye™ device and Ribo-Ker™ solution are CE marked.

The C-Eye™ device and Ribo-Ker™ solution are not available for sale in the United States.

#### Distributor:

Mediconsult AG  
Frohheimstrasse 2  
9325 Roggwil TG  
Schweiz  
071 454 70 20

Mediconsult SA  
Route du Pâqui 1  
1720 Corminboeuf  
Suisse  
026 467 60 20



[info@mediconsult.ch](mailto:info@mediconsult.ch)  
[www.mediconsult.ch](http://www.mediconsult.ch)

# ONE DEVICE TO PERFORM BOTH...



For a complete list of compatible slit lamps and distributors, visit [www.emagine-eye.com](http://www.emagine-eye.com)



## CXL at the slit lamp

The modern way



## CXL in the operating room

The classic way - with a table mount



### Innovative

The C-Eye™ device brings cross-linking technology to the slit lamp. Simple, efficient and safe from the novice to the expert user.



### Versatile

Perform CXL at the slit lamp. Or use the C-Eye device just like any other CXL device: on a table mount, in your operating room.



### Portable

Work in different places? Take the C-Eye™ device wherever you go. Charge it via USB-C, or even use a powerbank on-the-go.

# RIBO-KER™

ONE RIBOFLAVIN FOR KERATOCONUS AND KERATITIS TREATMENT

EMAGine is proud to offer Ribo-Ker™, its exclusive CE-marked riboflavin.

No more confusion with different types of riboflavin: EMAGINE has developed one solution for all applications.

- ✓ 0.1% hypoosmolaric riboflavin
- ✓ CE marked
- ✓ Ready for all epi-off CXL and all PACK-CXL applications
- ✓ Can be used on normal thickness and thin corneas
- ✓ Comes in a procedure pack with a speculum and a consumable cap for the C-Eye device



[www.emagine-eye.com](http://www.emagine-eye.com)

The C-Eye™ device and Ribo-Ker™ solution are CE marked.  
The C-Eye™ device and Ribo-Ker™ solution are not available for sale in the United States.

**Distributor:**  
Mediconsult AG  
Frohheimstrasse 2  
9325 Roggwil TG  
Schweiz  
071 454 70 20

Mediconsult SA  
Route du Pâqui 1  
1720 Corminboeuf  
Suisse  
026 467 60 20

mediconsult

[info@mediconsult.ch](mailto:info@mediconsult.ch)  
[www.mediconsult.ch](http://www.mediconsult.ch)