

OCULUS GONIO ready® CL
Self-Stabilizing, Direct
Gonioscopy Lens



Visualization of the Anterior Chamber Angle during MIGS


Innovative Design Ensures Excellent Lens Stability

The GONIO ready® CL offers the surgeon a crystal-clear view of the anterior chamber angle during MIGS without having to hold a lens. Due to its compact size and highly absorbent foam fin, this self-stabilizing gonioscopy lens ensures space for incisions while maximizing stability. As a single-use tool, the GONIO ready® CL is a timesaving, safe, and reliable addition to the OR.

Lens on, Hands off!


- Ultra high clarity through OCULUS quality lens
- Less adjustment due to large self-stabilizing foam fin
- Time saving: fast setup and easy use
- Improved safety: stable and sterile
- Ergonomic design for both left- and right-handed use

Great Stability of the Lens Due to Large Contact Surface



*Pictures: Prof. Dr. med. Marc Schargus, MHBA, FEBO
Head of Ophthalmology Department at Asklepios Klinikum Nord, Hamburg, Germany*

Brilliant View Through OCULUS Quality Lens



Technical Data

Article No.	GONIO ready® CL	box including 10 pieces
10012720		



WWW.OCULUSSURGICAL.COM



The OCULUS QM system is certified in accordance with ISO 13485 (MDSAP) and (EU) 2017/745 (MDR)

OCULUS Optikgeräte GmbH
Postfach • 35549 Wetzlar • GERMANY
Tel. +49 641 2005-0 • Fax +49 641 2005-299
Email: support@oculussurgical.de • www.oculussurgical.com

Find your local OCULUS representative on our website.

