

Dermapen 4™



One Device.
Tailored Solutions.
Endless Possibilities.

WORLD'S #1
MICRONEEDLING
COMPANY

Dermapenworld™

Upgrade to the Gold Standard in Microneedling: Dermapen 4™

More than a decade ago, DermapenWorld™ brought the original Dermapen™ into clinics. It revolutionised the microneedling industry as the world's first automated microneedling device.

Since then, our network of providers has delivered over 5 million Dermapen Treatments™ worldwide. DermapenWorld™ is now recognised as the World's #1 Microneedling Company.



WORLD'S #1
MICRONEEDLING
COMPANY

Dermapenworld™

Get the All-in-one Business Model

The Dermapen 4™ is a sleek, easy-to-use device that meets the highest design and production standards.

DermapenWorld supports serious practitioners who want to stay ahead. With Dermapen 4 in your clinic, you can offer advanced microneedling for any skin type, colour or concern with confidence, service, and outstanding patient results.

Is it the right choice for you? Discover why leading clinics choose the Dermapen 4.



Reasons To Upgrade to Dermapen 4™

Precision



Scar Treatment Setting to dramatically improve the look of scars.



AOVN™ Retention Valves for exact needle depths and reduced drag.



Built-in RFID Needle Cartridge ensures precise, reliable and drag-free needle penetration.



Sleek Digital Display to view and change procedural settings quickly and easily.

Quality



16 Needle Cartridge Design creates more rejuvenating channels for increased results with more comfort.



Dual Power Operation for ultimate convenience without compromising capability.



Superior AOVN™ 4 Technology for increased speed, comfort, and safety.

Safety



Anti-Contamination Management ACM™ for unrivalled safety and protection guaranteed.



Ergonomic Design means effortless, smooth, and comfortable procedures.

Business Potential



Offer Tailored Solutions with One Device by targeting multiple skin concerns at once.



Bluetooth Enabled Connectivity to keep your device updated and running smoothly.



Visionary Marketing Support to promote your practice and help attract more patients.



Secures your long-term investment.



Scar Treatment Setting

Dramatically Improve the Look of Scars

Scar revision represents a lucrative market you don't want to miss out on.

Dermapen 4 is the only microneedling device with a dedicated Scar Treatment (ST) setting and protocols for the revision of deep, atrophic scars, including surgical, burn contracture, striae and post-acne scars.

Needle penetration can reach up to 3.0mm to effortlessly and rapidly treat body contours and thickened, fibrous tissue. Needle depth may also be adjusted in 0.1mm increments for scalp, eyelids and intimate body areas.

Advanced Oscillation vertical Needle™ (AOVN) technology automatically calibrates the motor speed, pressure and penetration for targeted precision. There's no dragging. No pulling of the skin.

- Dermapen 4 is the only microneedling device with a dedicated Scar Treatment (ST) setting.



Digital Display. One Touch Workstation

View & Change Procedural Settings Easily

The illuminated LED display shows procedural depth, speed and battery status. You can change needle depth and treatment parameters with the simple push of a button.

Custom Dermapen 4 sterile sleeves further support a sterile working environment without compromising screen quality and clarity. During each procedure, the Digital-Servo motor accurately fine-tunes treatment parameters (including depth, speed and setting) to deliver exceptional accuracy.

Devices that use manual adjustments can't match the precision of the Dermapen 4, which can fine-tune needle depth in tiny 0.1mm steps. Plus, its built-in auto-calibration ensures consistent accuracy, eliminating any variations during treatments.




Built-in RFID Needle Cartridge

Ensures Effective & Exact Procedures

This patented breakthrough RFID technology guarantees that the desired depth is exact and ensures precise, reliable and drag-free needle penetration.

By automatically calibrating the needle cartridge with the device, all variations in production or wear and tear are accounted for. This ensures the reliability and consistency of treatments for the lifetime of your device.



AOVN™ Retention Valves

Exact Needle Depths & Reduced Drag

Specialised retention valves relieve the build-up of procedural pressure within the Dermapen 4 needle cartridge. The effect prevents fluid build-up within the needle cavity.

The vacuum-like sensation removes the risk of cross-contamination while ensuring pressure-free contact to the skin during a procedure.

As a true depth regulator, the AOVN vacuum-retention valves promote procedural accuracy. This means the depth displayed on the digital screen is the true measure of needle penetration into the patient's skin.

Regulated suction release ensures up to 1,920 micro-channels are smoothly delivered every second with no drag or resistance.

Needle sharpness and integrity are never compromised, even on large treatment areas.



Revolutionary 16 Needle Cartridge

Faster Treatment Times.
Larger Procedural Areas.

Generating 1,920 puncture channels per second, Dermapen 4's visionary 16-needle cartridge generates over 100% more puncture channels than counterfeit needling pens for faster treatment times and enhanced patient comfort.

The 33 gauge, medical-grade stainless steel needles are gamma-sterilised and oscillate up to 120 revolutions per second, equalling an unrivalled 7,200rpm. Only Dermapen 4 initiates more pathways to facilitate skin regeneration and the infusion of active substances. No other device on the market today comes even remotely close.

With patented technology, the Dermapen 4 needle cartridge offers a splayed and scalloped surface that glides over the skin. When combined with Dp Dermaceuticals Meso-Glides™, it offers a superior infusion of active ingredients. Testing has shown these micro-channels carry as much as 80% more topical nutrients into the skin.

The ISOgeometric placement of all 16 needles within a single cartridge head was meticulously designed by a team of industrial designers. After months of rigorous testing, this optimal arrangement ensures the needles glide smoothly over the skin, delivering exceptional results.



Dual Power Operation

Ultimate Convenience Without
Compromising Capability

Dermapen 4 offers an Advanced Lithium Turbo-Cut battery system that performs at the same level as mains power.

Dermapen 4 can provide up to 90 minutes of true operating time at full power capacity with just one battery. One battery also provides over 8 hours of standby.

Unlike traditional batteries, which gradually lose their power, causing uneven results and dragging of the skin, the Advanced Lithium Turbo-Cut technology ensures full and consistent power for the life of the battery.

Dermapen 4 also operates on traditional mains power with a 3.5-metre cord for added convenience.



Advanced Oscillating Vertical Needle™ (AOVN) Technology

Unbelievably Fast & Powerful Microneedling

The Dermapen 4™ needles puncture the skin in an up-and-down motion for increased precision, speed and safety.

MORE SPEED
A direct digital drive outperforms spring-loaded technology. Faster needles = faster procedures.

MORE POWER
Automated coupling of the cartridge to the Dermapen 4 handpiece maximises power.

MORE PRECISION
Digitally controlled needle movement enhances your accuracy for better results.

MORE SAFETY
The oscillation mechanism contains both internal and external fluid seals.

THE FINAL RESULT?
More impressive results and satisfied patients who refer others to your practice.

Dermapen 4's visionary
16-needle cartridge
generates over 100% more
puncture channels than
counterfeit needling pens for
faster treatment times and
enhanced patient comfort:





Ergonomic Design

Treat Face, Body, Scalp with Ease,
Comfort & Safety

The sleek design of Dermapen 4 creates an ergonomic weight distribution. Fitting easily into the hand, less pressure and grip are required to achieve needle contact on the skin.

The comfortable position removes practitioner fatigue and allows for multiple consecutive procedures. With less pressure, patients experience increased procedural comfort whilst maximising their results.

Dermapen 4 is so light and ergonomic it fits like a weightless glove.



Anti-Contamination Management™ (ACM™)

Offer the Safest Device Ever

Only Dermapen 4™ gamma-sterilised needle cartridges feature a unique, medical-grade silicone seal and fluid guard. The external cartridge is fortified with a specialised fluid guard that prevents cross-contamination to the outer cartridge shell.

This expertly engineered, patented, anti-contamination technology prevents the internal and external backflow of procedural fluids, even during the most challenging treatments.

ACM technology is a world-first and is exclusive only to Dermapen 4. Guaranteeing no cross-contamination risks for optimised patient safety.

The comfortable position
removes practitioner fatigue
and allows for multiple
consecutive procedures:





Offer Tailored Solutions with One Device

Target Multiple Skin Concerns at Once

Dermapen 4™ microneedling promotes rapid lifting, firming and toning of facial features and contours. It evens out skin tone and texture and diminishes scars and stretch marks. The overall result is a smoother, younger-looking complexion.

You can target multiple skin concerns at once to help any patient with confidence—no matter what skin type, colour or concern they come to you with.

It's your all-in-one business model for any skin type or concern.

PROFESSIONALLY ENDORSED

Thousands of doctors endorse our world-leading technology that helps your patients look younger quickly, without side effects or downtime.

PROVEN PROTOCOLS

Exclusively available to our Authorised Treatment Providers, these protocols are step-by-step guides to building a true Dermapen Treatment™.

SYNERGY OF SOLUTIONS

Get everything your clinic needs to thrive in one place. All our procedures and products are designed to work together for optimum results.



Marketing Support

Attract More Patients into Your Clinic

BE RECOGNISED

Receive an Authorised Treatment Provider plaque, exclusive marketing materials, point-of-sale tools, and clinic certification. This recognition helps build your clinic's credibility and attract more clients.

PRIORITY LISTING

Only Dermapen 4 providers get priority listing on DermapenWorld's Clinic Finder. Patients looking for genuine Dermapen Treatments will be directed to your door.

FREE EXCLUSIVE MEMBERSHIP

When you register your device, you'll also gain access to the dedicated Professional Portal. Discover a world of step-by-step Protocols and training, as well as marketing support to help you get the best commercial results.



Bluetooth Enabled Connectivity

Automatic Updates, Never Becomes Obsolete

Your device stays current because the software will be updated for you as new technologies are unearthed. The DermapenWorld App is easy to use, featuring:

- Completely reworked app architecture
- Use one app for Dermapen 4™ and Dermapen HOME™
- App screens and functions tailored to your account and requirements.
- Step-by-step guided set up for new users
- Instructional videos on each set-up page.



7 Year Warranty

Secures Your Long-term Investment

Including parts and labour, a 7-year, fully transferable and unrestricted warranty supports Dermapen 4 - the world's most trusted - "gold standard" device.

- NO OTHER PEN offers faster or more superior treatment results.
- NO OTHER PEN allows you to tap into the scar treatment market and maximise your ROI.
- NO OTHER PEN connects practitioners with a worldwide network of support and education.
- NO OTHER PEN competes with the safety and protection offered by Dermapen 4.
- NO OTHER PEN offers a 7-year, unconditional warranty.

Target multiple skin concerns
at once to help any patient
with confidence—no matter
what skin type, colour or
concern they come to you with.

Dp DOCK™ Compatible

Designed with skincare professionals in mind, the new Dp DOCK™ offers a Sync & Go concept that simplifies Dermapen Treatments™. This innovative docking station ensures your Dermapen 4™ device is always up to date with the latest firmware, allowing you to focus on delivering exceptional care to your clients.

The Dp DOCK™ features clear battery level indicators, convenient device storage, and integrated battery charger connections, allowing you to monitor your device at a glance. This ensures your tools are always ready for use, enhancing your operational efficiency and patient satisfaction.



Upgrade to Unlimited Possibilities with Dermapen 4™

Experience the fastest, most powerful microneedling device on the market. Thanks to superior precision, quality and safety features, treatments are completed in record time, and discomfort is minimised for patients.

But that's just the beginning...

Because it's suitable for all skin types, colours and concerns, the possibilities are endless. You can help any patient reach their skin goals faster with confidence thanks to Dermapen 4™

Say goodbye to complexity with the ultimate glow-up solution. Dermapen 4™ targets wrinkles, fine lines, acne, pigmentation, rosacea, and more—simultaneously.

**Are you ready to upgrade to Dermapen 4™
from the World's #1 Microneedling Company?**

“ The Dermapen 4 microneedling device offers unrivalled precision, minimising tissue trauma and allowing me to safely treat close to the eye. With a speed of 1920 micro-needles per second, we can treat a full face in as little as 10 minutes with minimal discomfort. Our patients benefit from virtually pain-free sessions without the need for numbing cream, minimal downtime and outstanding results for skin rejuvenation, pigmentation and scar remodelling, even after just one treatment. ”

Dr. Rachna Murthy ★★★★★
Consultant, Oculoplastic & Aesthetic Surgeon

WORLD'S #1
MICRONEEDLING
COMPANY

Dermapenworld™





Dermapenworld

Distribution in Switzerland:



LASERMED SA - Innovating Medicine

info@lasermed.ch | www.lasermed.ch

Frohheimstrasse 2 | 9325 Roggwil TG Rte André-Piller 18 | 1762 Givisiez FR
Tel. 071 454 70 30 | Fax 071 454 70 34 Tél. 026 466 38 15 | Fax 026 466 38 16